

**Small Anthropometer  
Model 01291**

[info@si-instruments.com](mailto:info@si-instruments.com)  
[www.si-instruments.com](http://www.si-instruments.com)



**Product Description**

Our Small Anthropometer has a range of 0 to 30 cm in 0.1 cm increments. Popular uses include measuring wrist, elbow, knee and ankle widths, as well as measuring smaller muscle masses like the bicep and calf. Aluminum in construction, it uses a spring-loaded ball bearing in a sliding C-shaped arm to provide accurate and precise measurement.