

## Two-Point Aesthesiometer Model 16011



### Product Description

Useful in determining the two-point threshold of touch. Easy to read slide scales are calibrated to the nearest 0.1 cm. Special vinyl tips minimize influence of temperature on sensation.

### Related Products

#### Three-Point Aesthesiometer

Model 16012

The third point allows easy alternation between one and two-point testing.

---

### Contact Information

S.I. Instruments  
256 South Rd. Hilton  
South Australia 5033  
Ph (08) 8352 5511

[info@si-instruments.com](mailto:info@si-instruments.com)  
[www.si-instruments.com](http://www.si-instruments.com)